

# SUPINE STABILISER

## Instructions for Use:

Place the stabiliser so the widest part is on the bottom, resting on the sheet. Support the person's thighs and place their legs through the gaps in the middle. The thin part of the pommel must be nearest the person's torso. Ensure the legs are fully supported.

Our Supine Stabiliser is an innovative new piece of equipment which helps to stabilise the femora and pelvis when a person is in a supine position. It is made from a breathable, washable fibre which draws vapour and sweat away from the body. This allows the individual using the Stabiliser to sleep in a temperate and hygienic environment

Note: 4 sizes available:


size 1 (infant), size 2 (18mths - 3 yrs), size 3 - (child to 14 yrs) size 4 (teenager - adult)



**C1SOUTH**

0800 212 333 email: [info@c1south.co.nz](mailto:info@c1south.co.nz) [www.c1south.co.nz](http://www.c1south.co.nz)

# SIZING GUIDE



**Size 1**

A	340mm
B	20mm
C	100mm
D	60mm
E	80mm

**Size 3**

A	420mm
B	50mm
C	180mm
D	100mm
E	90mm

**Size 2**

A	340mm
B	40mm
C	160mm
D	70mm
E	100mm

**Size 4**

A	480mm
B	60mm
C	220mm
D	130mm
E	110mm