



GLIDE WEAR^{by Tamarack}



Dual-layer Shear Reduction
Zone glides with users skin



Single layer stability zone enhances
user positioning and control

GLIDEWEAR CUSHION COVER


- Glidewear Seating Interface cushion cover slips over most wheelchair cushions
- Features dual-layer fabric technology that slides against itself for extra protection from friction and shear
- Shear Reduction Zone absorbs friction and shear forces under bony prominences
- Helps maintain proper positioning and stability while sitting
- Lightweight, breathable fabric helps reduce moisture build-up to preserve skin integrity

3 SIZES

Small 14x14 to 16x16

Medium 16x16 to 18x18

Large 18x18 to 20x20



*glide*free
Transfer Board Strips with ShearBan® Technology

GLIDEFREE TRANSFER BOARD TAPE KIT

- Provides a durable and low-friction surface
- Can be applied to nearly any wood, plastic or composite materialised-based transfer board
- Reduces surface friction to make sliding transfers easier
- Protects sensitive skin from abrasion and skin breakdown
- Can be used for wet and dry transfers

C1 SOUTH

AVAILABLE TO BUY ON WWW.C1SOUTH.CO.NZ