


WHEEL BLADES

GLIDE THROUGH WINTER WITH THE AWARD WINNING **WHEELBLADES S** WHEELCHAIR SKIS



reddot design award
winner 2013



Scan with your smart
phone for a video demo

Wheelblades S are an award winning, high-end mobility ski that attach onto the small front wheels of a manual wheelchair, allowing wheelchair users to glide through ice and snow, and improve mobility during winter. Sold as a pair.

The wide contact surface evenly distributes pressure onto the ground and prevents the small front wheels from sinking into the snow.

Wheelblades S are simple and quick to use with an adjustable clamp lock to fit wheels 1.8cm to 6cm wide. The lightweight design is portable and easy to install.

Installation: Simply place the Wheelblades S on the ground, tilt the wheelchair slightly and insert the small front wheelchair wheels into the binding shaft. Lock the Wheelblades clamp fastener and you are ready to go.



WHEEL BLADES

GLIDE THROUGH WINTER WITH THE AWARD WINNING **WHEELBLADES XL** STROLLER SKIS



XL



Scan with your smart
phone for a video demo

Wheelblades XL are an award winning, high-end ski that attach onto the front wheel(s) of a baby stroller for increased mobility through ice and snow. Sold as a pair or a single unit to accommodate 3 and 4 wheel strollers.

The wide contact surface evenly distributes pressure onto the ground and prevents the front wheels from sinking into the snow.

Wheelblades XL are simple and quick to use with an adjustable clamp lock to fit wheels 2cm to 7cm wide. The lightweight design is portable and easy to install.

Installation: Simply place the Wheelblades XL on the ground, tilt the stroller slightly and insert the front wheels into the binding. Lock the clamp fastener and you are ready to go.



C1SOUTH

0800 212 333

email: info@c1south.co.nz

www.c1south.co.nz

Jul 2015