


Dual-layer Shear Reduction Zone™
glides with users skin.



Single layer stability zone enhances
user positioning and control

The GlideWear Seating Interface cushion cover slips over most wheelchair cushions, providing extra protection from friction and shear.

- **Shear Reduction Zone:** Dual-layer fabric technology designed to slide against itself, absorbing friction and shear forces under the bony prominences of your pelvis. Increases the margin of safety where skin is most susceptible to breakdown.
- **Positioning and Stabilisation:** Friction is retained in the areas surrounding the 'Shear Reduction Zone' with a unique design within the interface to help maintain proper positioning and stability while sitting.
- **Versatile:** Designed to easily wrap around most wheelchair cushions.
- **Breathable:** Lightweight, breathable fabric helps reduce moisture build-up to preserve skin integrity.
- **Construction:** 85% Nylon. 15% Spandex. 100% Polyester Thread.

Product Code	Cover Size	Compatible Cushion Sizes
GW-SIT-SML	Small	14x14 to 16x16
GW-SIT-MED	Medium	16x16 to 18x18
GW-SIT-LRG	Large	18x18 to 20x20