



Mattress Overlay

For the prevention and treatment of pressure ulcers





Phone: 0800 212 333

 Antimicrobial agent incorporated which is effective against a wide range of micro-organisms including MRSA and E.coli.





Mattress Overlay

Clinically proven static air technology

Redistribution of pressure contours and supports around bony prominences while the ventilation holes allow for air flow which results in comfort and pressure injury prevention and treatment.

For home and hospital style beds






Sizing Guide (When Inflated)

Code	Size (cm)	Weight Limit
H5 4000A	193 × 86 × 8 (single bed size)	270 kg
H5 4010A	193 × 120 × 8 (bariatric)	360 kg





All come standard
with Double Action
Air Delivery Pump
for accurate
inflation.