




"Leaving a Lasting Impression"

No Weight Limit



Standard Sizing Options

Cushion Size	Code for 4" cell height	Code for 2" cell height	# of cells
14 x 14	SL4-1313-1	NA	7 x 7
14 x 16	SL4-1315-1	NA	7 x 8
14 x 18	SL4-1317-1	NA	7 x 9
16 x 16	SL4-1515-1	SL2-1515-1	8 x 8
16 x 18	SL4-1517-1	SL2-1517-1	8 x 9
16 x 20	SL4-1519-1	SL2-1519-1	8 x 10
18 x 16	SL4-1715-1	SL2-1715-1	9 x 8
18 x 18	SL4-1717-1	SL2-1717-1	9 x 9
18 x 20	SL4-1719-1	SL2-1719-1	9 x 10
20 x 16	SL4-1915-1	SL2-1915-1	10 x 8
20 x 18	SL4-1917-1	SL2-1917-1	10 x 9
20 x 20	SL4-1919-1	SL2-1919-1	10 x 10

The Starlock technology offers the ultimate in aggressive positioning, unparalleled stability with the highest level of therapeutic prevention and treatment of pressure ulcers and skin integrity issues.

Features:

- Available in 2" and 4" cell height.
- 'Lock' technology allows for easy adjustment to manage asymmetrical issues.
- 'Lock' technology allows you to isolate air cells and customise with "no contact" areas without the added expense of custom cushions.
- Ability to isolate airflow between cells enhances stability and assists with transfers and overall postural control.
- Maintains skin integrity by providing a low pressure, low shear and low friction seating environment.
- · Light weight approx 1.36kg.
- Comes complete with fitted cover, hand pump and repair kit.
 Alternative covers are available on request e.g. incontinent.
- Designed of high quality flame resistant latex free neoprene rubber.
- 24 Month manufacturers limited warranty.



Bariatric Sizing Options

Cushion Size	Cushion Code	# of cells
20 x 16	SL4-1915-1	10 x 8
20 x 18	SL4-1917-1	10 x 9
20 x 20	SL4-1919-1	10 x 10
22 x 18	SL4-2117-1	11 x 9
22 x 20	SL4-2119-1	11 x 10
24 x 18	SL4-2317-1	12 x 9
24 x 20	SL4-2319-1	12 x 10
26 x 20	SL4-2519-1	13 x 10
26 x 22	SL4-2521-1	13 x 11

C1SOUTH

Feb 2013

0800 212333 email: info@c1south.co.nz www.c1south.co.nz

***STARLOCK** CXR Cushion



CXR® Cushion Features:

- The CXR is designed for those individuals who require added protection and treatment under the coccyx area. Ideal for the individual who has, or is at risk of pressure ulcers.
- Comes standard with Starlock technology giving you the unprecedented stability and positioning that the Starlock cushion provides but with increased protection for the coccyx/sacral area.
- Light Weight, Approximately 3 pounds/1.36 kilograms.
- 24 Month Manufacturers Limited Warranty.
- Comes complete with fitted cover, hand pump and repair kit.
- Designed of high quality, flame resistant, latex free neoprene
- No Weight Limit!



Innovative Starlock technology creates a multitude of zones offering unprecedented positioning and stability. The Starlock set up is quite simple...adjust the Starlock like a traditional single valve air cushion. Then using the attached bulb pump, inflate the locking mechanism to secure client position. When locked the amount of air volume created in each individual air cell by virtue of the set up procedure will remain... creating a very personalised seating surface without the cost of a custom cushion.

Skin Protection, Positioning and Stability



*Configuration based on a 54 celled air cushion. (8 x 8 array)

CXR Sizing Chart

Cushion Size	Code	# of cells
16 x 16	CXR-1515-1	8 x 8
16 x 18	CXR-1517-1	8 x 9
16 x 20	CXR-1519-1	8 x 10
18 x 16	CXR-1715-1	9 x 8
18 x 18	CXR-1717-1	9 x 9
18 x 20	CXR-1719-1	9 x 10
20 x 16	CXR-1915-1	10 x 8
20 x 18	CXR-1917-1	10 x 9
20 x 20	CXR-1919-1	10 x 10