




- Multi-layered contoured foam cushion utilising the latest Visco Gel foam technology.
- The Visco Gel beads act to reduce heat build up, critical to minimise skin integrity issues. Made with antibacterial properties.
- **Selected Perforated Softening (SPS)** at the pelvic area enhances immersion for bony prominences. Also improves pelvic positioning, enhances comfort and stability.
- Designed for positioning and comfort for users who are at moderate to high risk of skin breakdown.
- Lightweight >1kg
- Comes in sizes from 10" to 20" wide.
- Internal and external covers provided as standard.

A = 4.25" (3.5" for 12" and 14" wide)

B = 3.75" (3" for 12" and 14" wide)

C = 3.25" (2" for 12" and 14" wide)



Code	Size*	Weight Limit
NXCABF 1010	10 x 10	80kg
NXCABF1212	12x12	80kg
NXCABF1414	14x14	100kg
NXCABF1416	14x16	100kg
NXCABF1616	16x16	136kg
NXCABF1618	16x18	136kg
NXCABF1620	16x20	136kg
NXCABF1816	18x16	136kg
NXCABF1818	18x18	136kg
NXCABF1820	18x20	136kg
NXCABF2016	20x16	136kg
NXCABF2018	20x18	136kg
NXCABF2020	20x20	136kg
NXCABF2220	22x20	136kg

*Please note that the cushion depth is taken from the side of the cushion due to the radius cut at the back of the cushion to account for upholstery or back supports.



Anatomically shaped SPS (Selectively Perforated Softening) perforations are anatomically positioned to support the Ischial/sacral area to prevent sliding and shear and maintain pelvic positioning. Provides immersion and lower pressures to the bony prominences while supporting surrounding areas that can take a higher load.

